



NHS SCOTLAND BROCHURE

PATIENT FURNITURE

CONTRACT REFERENCE: NP785/23

About Renray

Renray Healthcare is an approved NHS supplier with over 50 years' experience of manufacturing and supplying high-quality furniture including hospital patient and visitor chairs, lockers, overbed tables and more to the NHS.

Many products have been designed to meet the requirements of modern healthcare and hospital environments, incorporating design features to address infection control issues, reduce manual handling and pressure sore risks, and enhance patient comfort and well-being.

The increasing need for hospital furniture to withstand the stresses and strains of everyday life is imperative within the NHS. With rigorous research and development, we are constantly meeting this demand with our ranges of furniture.

We manufacture and assemble our healthcare furniture in our own purpose built factory in Winsford, Cheshire, as well as in Europe to British Standards. Hence we are able to ensure that your furniture is produced to the highest standards and we work with you to arrange, plan and meet your time schedule and budget.

Whether you are looking for a complete revamp of your hospital ward or looking to add some colour to your waiting area, our experienced sales staff are available to offer comprehensive advice on the most suitable furniture and upholstery fabrics to meet your needs.

Contract Reference:
NHS Scotland Contract NP785/23



Contents

FABRICS & FINISHES	4 - 5
PATIENT OVERBED TABLES	6 - 9
PATIENT LOCKERS	10 - 11
SEATING	12 - 18
- Balero	12
- Finland	14
- Seating Options	15
- Shrewsbury	16
- Kingsley	17
- Foot Stools	18
TABLES	20-21
- Coffee Tables	20
- Dining Tables	21
CONTACTS	22 - 23

For more information or to discuss upcoming projects contact NHS Scotland's dedicated Regional Sales Manager:

Mark Rodger
Mark.Rodger@renrayhealthcare.com
07780 955 225

High performance fabrics & finishes for a ward environment

Amongst our extensive range of vinyl and impervious fabrics, we offer a number of anti-microbial options that are effective against MRSA as well as E.coli, Listeria, Salmonella, Legionella and other similar bacteria.

Vinyl Fabrics



Impervious Fabrics



Please request a sample fabric swatch or a copy of the cleaning & care instructions on your chosen fabrics from our Sales Office on 01606 593456 or e-mail sales@renrayhealthcare.com

Fabric Cleaning Guide

Source: http://www.panaz.com/media/site-downloads/CARE_MAINTENANCE_OF_UPHOLSTERY_updated_.pdf

Provided by Panaz Fabrics

Care and maintenance of vinyl upholstery

01. Wipe clean regularly with a microfibre cloth to remove dust particles and light soiling.
02. Remove spillages promptly with an absorbent dry cloth.
03. For general soiling, use a microfibre cloth (taking care to avoid abrasive damage to surface) and tepid soapy water (non-alkaline).
04. Rinse clean with water and dry with an absorbent cloth.
05. Bodily fluids should be removed promptly within 15 minutes with cold water and then cleaned as above.
06. Extreme soiling/infection control purposes use a 5% sodium hypochlorite bleach solution, Haz Tabs or Chlorclean.
07. Steam cleaning is an effective sanitizer as are Haz Tabs and Chlorclean but is not a good cleaner. A wipe with a cloth, preferably microfibre, and water is more effective for actual cleaning. Steam and heat from steam cleaners actually seal in oily residues from bodily fluids, etc. If steam cleaning is essential, use infrequently for emergency infection control purposes only and in accordance with the machine manufacturer's recommendations. As we are unable to comment on the variety of machines in the marketplace or on the operator's handling, care should be taken for the sake of both the upholstery fabric and the wooden frame and we are unable to guarantee these.
08. Do not machine wash or dry clean.
09. Do not use furniture polish, stain remover or solvents.
10. Remove indigo dye from denim and canvas garments immediately as they will produce permanent staining if not dealt with promptly.
11. Do not use a felt tip or ball-point pen to mark.

Wood Finish Options



Pippy Oak
PO



Beech
B



Odessa
Golden Oak
OGO

To make an enquiry or to purchase this product call 01606 593456 or email sales@renrayhealthcare.com

Patient Over Bed Tables

Enterprise

The Enterprise™ tables are lightweight, durable and easy to clean. For complete flexibility they have been designed for the patient to use both in bed and when seated in a chair.




Enterprise tilting overbed table/chair with plastic base cover
PT11KAMBC-Beech

Code	Description	Table Top
PT11KAMBC	Enterprise Tilting Over Bed/ Chair Table w/ Moulded Top & Plastic Base Cover Finish - Beech Optional Finish: PT11KCMBC - Grey	



Enterprise non-tilting overbed table/ chair with plastic base cover
PT10KAMBC-Beech

Code	Description	Table Top
PT10KAMBC	Enterprise Non-Tilting Over Bed/ Chair Table w/ Moulded Top & Plastic Base Cover Finish - Beech Optional Finish: PT10KCMBC - Grey	

- 50mm twin wheel castors for easy manoeuvrability
- Gas spring assisted height adjustment
- Excellent height range of 710 - 1080mm
- Tilting top adjusted to a variety of angles using lightweight, easy to use ratchet system (PT11)
- Smooth, non-absorbent surfaces for ease of cleaning
- Table top 900mm x 420mm
- Safe working load (evenly distributed) 25kg
- Product weight 17kg
- Dimensions H710-1080 W900 D420mm

To make an enquiry or to purchase this product call 01606 593456 or email sales@renrayhealthcare.com

Patient Over Bed Tables

Elite

The Elite™ overbed table is lightweight, durable and easy to clean. For complete flexibility they have been designed for the patient to use both in bed and when seated in a chair.



- Smooth, non-absorbent surfaces for ease of cleaning
- 50mm twin wheel castors for easy maneuverability
- Gas spring assisted height adjustment
- Excellent height range of 710 - 1080mm
- Table top 900mm x 420mm
- Safe working load (evenly distributed) 25kg
- Product weight 17kg
- Dimensions H 710-1080 W900 D420


Code	Description	Table Top
E11V	Elite Tilting Overbed/ Chair Table Vinyl Wrap Top Finish: Pippy Oak	

Wolverhampton

The Wolverhampton provides a stable and smooth surface to eat, work or enjoy activities from.



- Fully moulded, smooth, non-absorbent surface for ease of cleaning
- Height adjustable table for use over your bed or armchair
- Minimal dirt traps
- Rounded corners
- 65mm castors for ease of mobility
- Dimensions H 680-1050 W380 D620

Code	Description	Table Top
WOLC	Wolverhampton Overbed Table Moulded Top Finish: Pippy Oak	

Patient Lockers

Type A4SC

Self medication lockers practically designed.



Layout:

- Lockable self medication drawer
- Middle drawer
- Cupboard door
- 50mm castors

Optional:

- 75mm braked castors
- Suited Locking System

Accessory:

MK Master Key Relacement

Code	Description	Height mm	Width mm	Depth mm
TYPE4SC	Patient bedside locker with self medicating drawer. Finish: Beech	850	550	530

Florida

The Florida range is without compromising on design, creating a sophisticated and contemporary look in its simplicity.



Features:

- 18mm MFC drawer and door fronts with PVC edging
- 18mm MFC tops with PVC edging
- Chrome bar handles
- Lockable bedside unit
- Available in a variety of wood finishes

Code	Description	Height mm	Width mm	Depth mm
FL116L	Florida Bedside Unit w/ Door & Lock Finish: Pippy Oak	660	445	457

To make an enquiry or to purchase this product call 01606 593456 or email sales@renrayhealthcare.com

Seating

Balero

Attractive and practical design with a robust frame



BL300 Balero High Back Chair
RS Removable Seat Cushion

- Sacral gap and an optional removable seat cushion allow easy cleaning and access
- An optional pressure reducing seat cushion for the management and prevention of pressure ulcers
- Optional housekeeping wheels make the chair easy to move reducing moving and handling risks
- All wooden frames are coated in lacquer and anti-microbial fabrics are available.

Go to page 15
for chair options

BL140 Balero Medium Back Chair
RS Removable Seat Cushion



BL105 Balero Low Back Chair
RS Removable Seat Cushion

Code	Description	Height mm	Width mm	Depth mm	Seat Height mm	Seat Width mm	Seat Depth mm
BL300	Balero High Back chair	1219	610	686	483 (19")	483	483
BL140	Balero Medium Back Chair	1041	610	635	483 (19")	483	483
BL105	Balero Low Back Chair	940	610	635	483 (19")	483	483

Optional Seat heights: 17", 18", 20"

Safe Working Load 25 stone

Petite Option: BL140+16, Seat Height 16"

To make an enquiry or to purchase this product call 01606 593456 or email sales@renrayhealthcare.com

Seating

Finland

A modern style available in a choice of contemporary colours.



Finland Stacking Chair w/ Arms
F835

- Sacral gap on all chairs for easy cleaning
- Strong and durable laminate frame
- Stackable 6 high
- Cushioned seat and back rest

F833 Finland Stacking Chair

Code	Description	Height mm	Width mm	Depth mm	Seat Height mm
F833	Finland Stacking Chair	850	550	530	19"
F835	Finland Stacking Chair w/ Arms	850	550	530	19"

Seating Options



House Keeping Wheels
H



Height Adjusters
J



Height Adjusters With Rear Wheels
JW



Drop Down Arms
DDA



Skids
S



Removable seat
Standard
RS



Removable Pressure Relief
Low/Med Risk Comfort
Cushion: VRS
Removable Pressure Relief
High Risk Pressure Relief
Cushion: VRSX

To make an enquiry or to purchase this product call 01606 593456 or email sales@renrayhealthcare.com

Seating

Shrewsbury

Shrewsbury chairs have been cleverly created to look good in any location.



Shrewsbury Tub Chair
SHRT

- Excellent durability and comfort
- Available in a variety of anti-microbial and waterproof upholstery for ease of cleaning
- Contemporary twist on a classic design
- Optional seat height 18"

Code	Description	Height mm	Width mm	Depth mm	Seat Height mm
SHRT	Shrewsbury Tub Chair	780	590	500	457

Kingsley

Affordable, contemporary & comfortable compact seating suitable for any living space.



Kingsley 2 Seater Sofa
KINL2S

- Suitable for traditional & modern living spaces
- Excellent durability and comfort
- Available in a variety of anti-microbial and waterproof upholstery for ease of cleaning
- Available as an armchair

Code	Description	Height mm	Width mm	Depth mm	Seat Height mm	Seat Width mm	Seat Depth mm
KINL2S	Kingsley 2 Seater Sofa	838	1219	711	457	1020	485

To make an enquiry or to purchase this product call 01606 593456 or email sales@renrayhealthcare.com

Seating

Foot Stools

Designed to complement the Renray patient seating range



FRO Royal Footstool



FAJ Height Adjustable Footstool

Code	Description	Height mm	Width mm	Depth mm
FRO	Royal Foot Stool	432	432	370
FAJ	Height Adjustable Foot Stool	369 - 494	510	305

Dining

Finland

A modern style available in a choice of contemporary colours.



F833 Finland Stacking Chair



Finland Stacking Chair w/ Arms
F835

- Sacral gap on all chairs for easy cleaning
- Strong and durable laminate frame
- Stackable 6 high
- Cushioned seat and back rest

Code	Description	Height mm	Width mm	Depth mm	Seat Height mm
F833	Finland Stacking Chair	850	550	530	19"
F835	Finland Stacking Chair w/ Arms	850	550	530	19"

To make an enquiry or to purchase this product call 01606 593456 or email sales@renrayhealthcare.com

Tables

Dining Tables

A variety of stylish shapes and designs to cover all dining requirements, in a smooth, heat resistant laminate finish.



DTC Circular Dining Table



DTS Square Dining Table



DTR Rectangular Dining Table

Code	Description	Height mm	Width mm	Depth mm
DTC	Circular Dining Table	760	1070 Dia	n/a
DTS	Square Dining Table	760	910	910
DTRS	Rectangular Dining Table 1200mm	760	800	1200
DTR	Rectangular Dining Table 1500mm	760	900	1500

* 600, 900, 1200 & 1500mm size options also available, please enquire

Coffee Tables

Versatile and stylish, available in an attractive smooth laminate finish.



CTC Circular Coffee Table



CTU Low Square Coffee Table

Code	Description	Height mm	Width mm	Depth mm
CTC	Circular Coffee Table	640	Dia 600	N/A
CTU	Square Coffee Table - Low Height	520	600	600

To make an enquiry or to purchase this product call 01606 593456 or email sales@renrayhealthcare.com

Contact Us

For further information, including product availability, the latest colour options and prices, to order any of the products in this brochure; or if you have any queries about other Renray products and services, contact:

Mark Rodger
Scotland Area Sales Manager
Mark.Rodger@renrayhealthcare.com
M: 07780 955225

John McEwen
Commercial Director
John.McEwen@renrayhealthcare.com
M: 07425 248 871

Alison Ballard
Internal Sales Contact
Alison.Ballard@renrayhealthcare.com
T: 01606 593456 Ext: 442

Christopher Taylor
Customer Service Executive
Christopher.Taylor@renrayhealthcare.com
T: 01606 593456 Ext: 433

Pam Roberts
Logistics Co-ordinator
Pam.Roberts@renrayhealthcare.com
T: 01606 593456 Ext: 426

Can't find what you're looking for?

Renray has a vast product portfolio so if you can't find an item to fulfill your requirement, view our wide range of furniture brochures by visiting renrayhealthcare.com/brochures

Alternatively, contact our
Sales Office on:
T. 01606 593456
sales@renrayhealthcare.com

NHS Scotland Contract References:
NHS Scotland Contract NP785/23



Road Five, Winsford Industrial Estate,
Winsford, Cheshire CW7 3RB

www.renrayhealthcare.com