

Patient and visitor furniture that makes you feel close to home.



WELCOME

About Renray

Renray Healthcare is an approved NHS supplier with over 50 years' experience of manufacturing and supplying high quality furniture including patient and visitor chairs, lockers, overbed tables and more to the NHS.

Many products have been designed to meet the requirements of modern healthcare and hospital environments, incorporating design features to address infection control issues, reduce manual handling and pressure sore risks, and enhance patient comfort and wellbeing.

The increasing need for hospital furniture to withstand the stresses and strains of everyday life is imperative within the NHS. With rigorous research and development, we are constantly meeting this demand with our ranges of furniture.

We manufacture and assemble our healthcare furniture in our own purpose built factory in Winsford, Cheshire, as well as in Europe, to British Standards. Hence we are able to ensure that your furniture is produced to the highest standards and we work with you to arrange, plan and meet your time schedule and budget.

Whether you are looking for a complete revamp of your hospital ward or looking to add some colour to your waiting area, our experienced sales staff are available to offer comprehensive advice on the most suitable furniture and upholstery fabrics to meet your needs.

NHS Contract References:

NHS Supply Chain Contract Reference: 2017/S 207-427709

Ward and Residential Furniture 2017/S 056-10455

Mental Health Furniture 2016/S 041-066764

NHS SBS Contract Reference SBS/19/MA/TYX/9363/37

NHS Scotland Contract NP785/19



Contents

	Patient bedside furniture		Seating		Beds, foams and PAC
06-07	Enterprise overbed & chair tables	40	Belton	62	Elite
08-09	Enterprise lockers	41	Forli	63	Prestige bariatric
10-11	Elite lockers	42	Kingsley	64	Carer
12-13	'A' Range lockers	43	Harvard	65	Corus
14	600 range lockers	44	Chelford	66	Mattresses
15	100 range combination lockers	45	Modica	66	- Corus mattress
16-17	Vision type A4 lockers	46	Teramo	67	- Nautilus mattress
18-19	Elite patient chair range	47	Milan	67	- Bariatric mattress
20-21	Prestige chair range			68	- Jet mattress
22	Balero chair			69	Cushions
23	National chair	50	Dining furniture	69	- Corus foam
24	Finland stacking chair	51	Dining tables	69	- Noodles technology
25	Pop plastic chair	52	Ridley table	69	- Aurora active
26	Titan antimicrobial chair	53	Chelford chair	70	PAC
27	Lamarr breastfeeding chair		Finland chair	70-71	- Noodles mattress
				72	- Dormir
				73	- Diamond
				74	Beds
				74	- Abbey divan bed
				75	- Denver box bed
				76	- Jetbed
				77	- Stretford divan bed
					Staff accommodation
				80	Florida range
				81	Albany range
					About Renray
				82	Fabrics and cleaning instructions
				83	Express delivery
				84	What we do
				85	Renray credentials
				86-87	Contact



Look out for the Express Delivery icon on many of our products.

See page 83 for more information.

Patient bedside furniture that makes you feel close to home.

Our high quality and practical patient furniture, designed for use in a ward setting, has been carefully constructed with unique options and to meet the needs of today's healthcare environments.



- 06-07 Enterprise tables
- 08-09 Enterprise lockers
- 10-11 Elite lockers
- 12-13 'A' Range lockers
 - 14 600 range lockers
 - 15 100 range combination lockers
- 16-17 Vision type A4 lockers
- 18-19 Elite patient chair range
- 20-21 Prestige chair range
 - 22 Balero chair
 - 23 National chair
 - 24 Finland stacking chair
 - 25 Pop plastic chair
 - 26 Titan antimicrobial chair
 - 27 Lamarr breastfeeding chair




Enterprise over bed & chair tables

Features:

- Smooth, non-absorbent surfaces for ease of cleaning
- Gas spring assisted height adjustment
- 50mm twin wheel castors for easy manoeuvrability
- Excellent height range of 710 - 1080mm
- PT11 - Tilting top adjusted to a variety of angles using lightweight, easy to use ratchet system
- Table top 900mm x 420mm
- Safe working load (evenly distributed) 25kg
- Product weight 17kg.

Colour options:



Code	Colour	Description	Height range (mm)	Table top (mm)	Width (mm)	Depth (mm)
PT10KA	Beech	Enterprise non tilting overbed table/chair. Plastic kydex top	710 - 1080	900 X 420	420	470
PT10KB	Blue					
PT10KC	Grey					
PT10KAMBC	Beech	Enterprise non tilting overbed table/chair with plastic base cover	710 - 1080	900 X 420	420	470
PT10KBMC	Blue					
PT10KCMBC	Grey					



PT10KB Enterprise non tilting overbed table/chair with plastic base cover



Code	Colour	Description	Height range (mm)	Table top (mm)	Width (mm)	Depth (mm)
PT11KA	Beech	Enterprise tilting overbed table/chair. Plastic kydex top	710 - 1080	900 X 420	420	470
PT11KB	Blue					
PT11KC	Grey					
PT11KAMBC	Beech	Enterprise tilting overbed table/chair with plastic base cover	710 - 1080	900 X 420	420	470
PT11KBMC	Blue					
PT11KCMBC	Grey					



PT11KA Enterprise tilting overbed table/chair with plastic base cover

The Enterprise™ tables are lightweight, durable and easy to clean. For complete flexibility, they have been designed for the patient to use both in bed and when seated in a chair.



Plastic base cover option



Gas spring lever



Enterprise lockers

The Enterprise™ lockers are extremely strong and durable, as well as having a novel contemporary appearance.



Features:

- One piece plastic moulded drawers and smooth surfaces for easy cleaning
- Bottle holder (holds 3 x 1.5L bottles)
- Free 65mm castors (2) and kick stop castors (2)
- Dual sided bag hooks and retractable towel rail
- Lightweight, easy to manoeuvre with ergonomic easily accessible drawers
- Sturdy corner buffers
- Large bottom drawer with internal shelf.


Colour options:



Beech



Lockable drawers



Braked castors



Base of locker



PL21B Enterprise locker, 2 drawer c/w 1 lockable drawer



Code	Colour	Description	Height (mm)	Width (mm)	Depth (mm)
PL21B	Beech	Enterprise locker, 2 drawer c/w 1 lockable drawer	850	450	470



PL31B Enterprise locker, 3 drawer c/w 1 lockable drawer



Code	Colour	Description	Height (mm)	Width (mm)	Depth (mm)
PL31B	Beech	Enterprise locker, 3 drawer c/w 1 lockable drawer	850	450	470

Optional	L-CAB	Lock	GMK	Grand master key - suited locking system
	MK	Master key	DL	Digital lock (3780)
	SL	Suited locking system	RF	RFID lock
	SMK	Suited master key		

Elite lockers



E1L Locker with lockable compartment and drawer



E1SC Locker with lockable compartment and door

The Elite locker is fully moulded with rounded corners, lockable medication compartments and minimal dirt gaps helping reduce infection control issues by being easy to clean. Its solid carcass will endure even the most demanding hospital environments, and its patented design gives it a modern edge.

Layout:


- Lockable self medication compartment
- Open compartment
- Cupboard door or drawer
- 55mm castors

Specification:

- Kydex plastic

Optional:

- Drawer
- 75mm braked castors



Code	Colour	Description	Height (mm)	Width (mm)	Depth (mm)
E1L	Grey/Blue	Locker with lockable compartment and drawer	950	450	450
E1SC	Grey/Blue	Locker with lockable compartment and door	950	450	450

Optional	MK	Master key	MKSLS	Master key - suited locking system
	SL	Suited locking system	GMSLS	Grand master key - suited locking system

Features:

- Applied patented design
- Moulded thermo-plastic panels
 - easy to clean
- Minimal dirt traps
- Surface of the plastic sheet is more robust than typical laminated wood panels
- Rounded corners
- No exposed internal fittings or fixings
- Standard colours (Blue and Grey).



'A' range lockers

Self medication lockers practically designed to offer clever storage solutions.



TYPE A2 Patient bedside locker with self medicating drop down flap



TYPE A2 Patient bedside locker with self medicating drop down flap



TYPE A3 Patient combination locker with self medicating drop down flap




TYPE A4 Patient bedside locker with self medicating drawer



TYPE A4SC Patient bedside locker with self medicating drawer



TYPEA5SC Type A Range combination locker



Code	Colour	Description	Height (mm)	Width (mm)	Depth (mm)
TYPE A2	Beech	Patient bedside locker with self medicating drop down flap	975	460	515
TYPE A3	Beech	Patient combination locker with self medicating drop down flap	1400	715	515
TYPE A4	Beech	Patient bedside locker with self medicating drawer	975	460	515
TYPE A4SC	Beech	Patient bedside locker with self medicating drawer	850	462	424
TYPEA5SC	Beech	Type A range combination locker	1350	722	507

Optional	BC	SMK	SL
	Braked castors	Suited master key	Suited locking system
	MK	GMK	
	Master key	Grand master key	

Features:

- 75mm 360 degree swivel castors (A2, A4)
- 50mm 360 degree swivel castors (A4SC, A5SC)
- Smooth laminate finish (A2, A4)
- Smooth MFC finish (A4SC, A5SC)
- Towel rail to the rear (A2)
- Bowl holder (A2)
- Lock with two individual keys (All).

Layout:

- Self medicating lockable compartment, drop-down flap, heavy duty hinge (A2, A3)
- Deep lockable top drawer (A4, A5SC)
- Through doors, accessible from the left and right (A2, A3)
- Middle drawer (A4, A4SC, A5SC)
- Bottom cupboard door (A5SC)
- Wardrobe on left (A5SC)
- 2 door cupboard with adjustable shelf (All)
- Robe 255 wide with hanging rail (A3)
- Self medicating compartment internal dimensions. w: 416mm x d: 146mm x h: 220mm
- Self medicating drawer dimensions. w: 354mm x d: 220mm x h: 319mm.

600 range lockers



DB602 Patient bedside locker self medicating drop down flap

Flexible design, made to a robust standard with a smooth laminate finish including a moulded top.

Features:


- Kydex moulded top
- 80mm 360 degree swivel castors and buffers
- Post formed smooth laminate sides
- Smooth laminate finish
- Towel rail to the rear
- Bowl holder optional
- Lock with two individual keys.

Layout:

- Self medicating compartment internal dimensions w: 416mm x d: 146mm x h: 220mm, drop down flap
- Through doors, accessible from the left and right
- Middle drawer
- 2 door cupboard with adjustable shelf
- Robe 255mm wide with hanging rail (DB652).



DB652 Patient combination locker self medicating drop down flap



Code	Colour	Description	Height (mm)	Width (mm)	Depth (mm)
DB602	Beech	Patient bedside locker self medicating drop down flap	1050	455	455
DB652	Beech	Patient combination locker self medicating drop down flap	1400	710	515

Optional	BC	Braked castors	SMK	Suited master key
	MK	Master key	GMK	Grand master key
	SL	Suited locking system		

100 range combi lockers



DB103 Combination locker 1 drawer and 1 door

High quality and versatile, covering various needs with many built-in features.

Features:

- 75mm 360 degree swivel castors
- Smooth laminate finish
- Chrome handle
- Towel rail to the rear (bowl holder optional).

Layout:


- 1 drawer/1 door cupboard with adjustable shelf (DB103)
- 4 drawers (DB101)
- Robe 255mm wide with hanging rail (DB101, DB103).



Braked castor



DB101 Combination locker 4 drawer



Code	Colour	Description	Height (mm)	Width (mm)	Depth (mm)
DB101	Beech	Combination locker 4 drawer	1400	715	515
DB103	Beech	Combination locker 1 drawer and 1 door	1400	715	515

Optional	BC	Braked castors

Vision type A4 lockers

The vibrant Vision Type A4 locker range was specifically designed for use by the NHS to add stylish colour solutions to create positive hospital environments.

Layout:


- Deep lockable self medication top drawer
- Middle drawer
- Cupboard drawer
- 65mm castors
- Smooth finish laminate top
- Chrome handle
- Choice of colourful laminate finishes.



Code	Colour	Description	Height (mm)	Width (mm)	Depth (mm)
VISA4+YEL	Sunshine Yellow	Bedside locker	850	462	422
VISA4+RAS	Refreshing Raspberry				
VISA4+PUR	Dusky Damson				
VISA4+PIN	Rose Pink				
VISA4+BL	Sky Blue				
VISA4+GR	Fresh Lime				
VISA4+GREY	Grey				
VISA4+B	Beech				

Optional	L-CAB	Lock	DL	Digital lock (3780)
	MK	Master key	SMK	Suited master key
	SL	Suited locking system	GMK	Grand master key

Colour options:



Elite patient chair range

Standard:



Pressure reducing visco cushion (V)



Removable seat (RS)



Housekeeping wheels (H)

Optional:



Seat board standard (S25) for use with dynamic cushion



Height adjusters (J)



Armpads (A)




Integral wing option



E301 Medium chair without wings



E600 Elite standard chair with wings



Code	Description	Height (mm)	Width (mm)	Depth (mm)	Seat height (mm)	Seat width (mm)	Seat depth (mm)	Max. weight
E600	Standard chair with wings	1110-1210	620	720	460-560	500	460	127kg/20st
E300	Medium chair with wings	1160	620		460			127kg/20st
E601	Standard chair w/o wings	1110-1210	575	720	460-560	500	460	127kg/20st
E301	Medium chair w/o wings	1160	575		460			127kg/20st


Height options: (E600) 16", 17", 18", 19", 20", 21", 22" (E601) 16", 17", 18", 19", 20" (E300) 18"

Moving and handling:

- The housekeeping wheels allow the chair to be easily moved (when chair is empty)
- Sacral gap allows easy access to fit hoist slings
- The arms allow easy grip with rounded ends and a stable base to push up from
- Splayed legs ensure greater stability and easy hoist access
- SWL 127kg (20 stone).

Tissue viability:

- The visco-elastic seat cushion helps in the treatment and management of pressure ulcers with frequent monitoring and repositioning
- Impressive pressure reduction testing
- The 10 degree recline assists spinal alignment
- Option of a seat board (SB) to allow specialist cushions to be used.



Infection control:

The Elite chairs have been specifically designed to be easy to clean.

Features include:

- Saddle stitched seams to reduce the risk of fluid ingress
- Removable seat cushion for optimal cleaning access
- Sacral gap to allow easy cleaning and access
- Optional anti-microbial fabrics
- Anti-bacterial lacquer coating on the wooden frame.

The Elite™ range has been designed using anthropometric data which highlights 13 critical dimensions to achieve a good sitting posture. These 13 dimensions are included in each chair size that make up the Elite™ range. Sitting in the correct sized chair ensures the patient's comfort and wellbeing.



Prestige chair range



A versatile, modern and comfortable patient chair, suitable for varying user weights.



Features:

- Sturdy arms allow easy grip with rounded ends and a stable base to push up from
- PTG725, PTG735, PTG750 - Adjustable drop down arms to aid lateral transfer and ensure flexible positioning for patients
- Sacral gap allows easy access to fit hoist slings
- Wooden frame is coated with antibacterial lacquer
- Easy to clean
- Available in a variety of user weights; 25, 35 and 50 stones
- 18" standard height
- Available with or without arms
- Available in range of fabrics and wood finishes to suit any area
- Removable seat as standard.



Code	Description	Height (mm)	Width (mm)	Depth (mm)	Seat height (mm)	Seat width (mm)	Seat depth (mm)	Max. weight
PTG325	High back chair	1210	605	675	457 or 508	490	480	159kg/25st
PTG335	Large high back chair	1210	675	675	457 or 508	560	480	223kg/35st
PTG350	Extra large high back chair	1210	820	675	457 or 508	705	480	318kg/50st



Code	Description	Height (mm)	Width (mm)	Depth (mm)	Seat height (mm)	Seat width (mm)	Seat depth (mm)	Max. weight
PTG725	High back chair with drop down arms	1210	605	675	457 or 508	490	480	159kg/25st
PTG735	Large high back chair with drop down arms	1210	675	675	457 or 508	560	480	223kg/35st
PTG750	Extra large high back chair with drop down arms	1210	820	675	457 or 508	705	480	318kg/50st



Code	Description	Height (mm)	Width (mm)	Depth (mm)	Seat height (mm)	Seat width (mm)	Seat depth (mm)	Max. weight
PTG150	Extra large low back chair	1000	820	675	457 or 508	705	480	318kg/50st



Code	Description	Height (mm)	Width (mm)	Depth (mm)	Seat height (mm)	Seat width (mm)	Seat depth (mm)	Max. weight
PTG050	Extra large low back chair no arms	1000	820	675	470	527	476	318kg/50st

Optional	V	Pressure reducing visco cushion	H	House keeping wheels
Optional	J	Height adjusters	JW	Height adjusters with rear wheels



PTG335 Large high back chair



PTG735 Large high back chair with drop down arms



PTG150 Extra large low back chair



PTG050 Extra large low back chair no arms

Balero chair



BL300 High back chair

Attractive and practical design with a robust frame.

Features:

- Sacral gap and an optional removable seat cushion allows for easy cleaning and access
- An optional pressure reducing seat cushion for the management and prevention of pressure ulcers
- Optional housekeeping wheels make the chair easy to move reducing moving and handling risks
- All wooden frames are coated in antibacterial lacquer and anti-microbial fabrics are available.



Code	Description	Height (mm)	Width (mm)	Depth (mm)	Seat height (mm)	Seat width (mm)	Seat depth (mm)	Max. weight
BL300	High back chair	1219	610	686	406 - 508	483	483	158kg/25st
BL140	Medium back chair	1041		128kg/20st				
BL105	Low back chair	940	533	635	406 - 508	483	483	128kg/20st
BL100	Low back chair w/o arms	940						

Seat height options: 16", 17", 18", 19", 20" - (19" as standard)

Optional	V	Pressure reducing visco cushion	A	Armpads
	RS	Removable seat cushion	W	Wings
	H	House keeping wheels	J	Height Adjusters
	JW	Height adjusters with house keeping wheels		



BL105 Low back chair BL140 Medium back chair

National chair



NAT601 High back chair with E4000SC pressure reducing seat cushion

The National chair with removable pressure reducing foam seat cushion in a two way vapour permeable polyurethane coated cover, ensures 24 hour pressure area care for patients.

Features:

- Removable pressure reducing seat cushion
- Castellated pressure relieving foam
- PU breathable cover to reduce lesions
- Sturdy arms allow easy grip with rounded ends and a stable base to push up from
- Housekeeping wheels, as standard
- Sacral gap allows easy access to fit hoist slings
- Wooden frame is coated with antibacterial lacquer
- Easy to clean
- 19" standard height.

Pressure reducing seat cushion colour options: Grey PU, Plum PU and Blue PU.



NAT600 High back chair with integral wings and E4000SC pressure reducing seat cushion



Code	Description	Height (mm)	Width (mm)	Depth (mm)	Seat height (mm)	Seat width (mm)	Seat depth (mm)
NAT601	High back chair with E4000SC pressure reducing seat cushion	850	462	422	482	502	457
NAT600	High back chair with integral wings and E4000SC pressure reducing seat cushion						

Seat height options: +16 = 16", +17 = 17", +18 = 18" and +20 = 20"

Optional	JW	Height adjusters with house keeping wheels	A	Armpads
	J	Height Adjusters		

Finland stacking chair


A modern style available in a choice of contemporary colours.

Features:

- Sacral gap on all chairs for easy cleaning
- Strong and durable laminate frame
- Recommended stacking six chairs maximum.



F834 Stacking chair with curved arms F835 Stacking chair with straight arms



Code	Description	Height (mm)	Width (mm)	Depth (mm)	Seat (h) (mm)	Seat (w) (mm)	Seat (d) (mm)
F833	Stacking chair	850	550	530	483	500	570
F834	Stacking chair with curved arms						
F835	Stacking chair with straight arms						

F833 Stacking chair

POP plastic chair

Part of Renray's Vision Colourful Ward Furniture collection. Strong, lightweight and durable one piece air moulded visitor chair, which makes it easy to clean. Suitable for stacking eight chairs high.



Code	Description	Colour	Height (mm)	Width (mm)	Depth (mm)	Seat height (mm)	Max. weight
POP+YEL	Pop chair	Sulfur Yellow	790	500	514	460	159kg/25st
POP+GREY	Pop chair	Signal Grey					
POP+WHITE	Pop chair	Traffic White					
POP+GREEN	Pop chair	Yellow Green					
POP+MARINE	Pop chair	Marine Blue					
POP+LILAC	Pop chair	Blue Lilac					
POP+CREAM	Pop chair	Cream					
POP+BLACK	Pop chair	Black					



POP Plastic chair

Features:

- 100% recyclable
- Lightweight and easy to handle
- No screws, joints or stress points
- Ergonomic with good lumbar support
- Suitable for outdoor use
- UV resistant
- Made from reinforced polypropylene
- EN 1729:2006 Parts 1 & 2 ergonomics and strength compliant
- EN 15373:2007 Level 2".

Titan antimicrobial chair

The only easy-clean seating range that is antimicrobial to the core. Clinically proven to kill >99.9% of MRSA and E.Coli as well as inhibit fungal growth.

Features:

- Easy to clean
- Designed to have no areas to trap dirt; following NHS guidance
- Enabling quick, effective all over cleaning, it can be safely cleaned using: detergents, steam cleaning and Actichlor Plus
- Antimicrobial to the core - the patented silver based antimicrobial additive is clinically proven to kill >99.9% of bacteria including MRSA and E.Coli on contact
- The anti-tamper construction has rounded edges and a stable anti-tilt base to prevent injury
- Our seating is designed to British Standards for ergonomics, to provide a user with support.



Code	Colour	Description	Height (mm)	Width (mm)	Depth (mm)	Seat (h) (mm)	Max. weight
T110	Sky Blue	Titan chair	897	635	535	460	158kg/25st
T100	Blue	Titan chair with arms					

T110 Antimicrobial adult chair - plastic



T100 Antimicrobial adult armchair - plastic

Lamarr breastfeeding chair

High backed breastfeeding chair with wider shoulders for more comfort and reduced lumbar width, to allow for baby in all feeding positions (cradle, cross-cradle, football etc.).



Code	Description	Height (mm)	Width (mm)	Depth (mm)	Seat (h) (mm)	Seat (w) (mm)	Seat (d) (mm)	Max. weight
LAMW	Lamarr breast feeding chair	970	780	740	480	520	560	254kg/40st



Features:


- Design of seat prevents mothers from postural pains and reduces pressure on the perineum
- Supports the mother and child in all breastfeeding positions without slipping
- Provides excellent support for arms, back and shoulders, whilst allowing for flexible manoeuvre of the child.



Waiting area furniture that makes you feel close to home.

These items are designed to complement the rest of our range to ensure that the overall look is consistent and there are appropriate pieces to fill every need. All have been made by our craftsmen with both style and durability in mind.

- 30-31 Waiting area coffee tables
- 32 Venetian chairs
- 33 Finland chairs
- 34 Waiting area tub chairs
 - 35 - Chelford
 - 35 - Kingsley
 - 36 - Balero
 - 36 - Teramo
 - 37 - Napoli
 - 37 - Shrewsbury



Waiting area coffee tables



Versatile and stylish, available in an attractive smooth finish with laquered edges.



CTS Square coffee table



CTU Square coffee table - low height



CTR Rectangular coffee table



CTC Circular coffee table



CTCL Circular coffee table - low height

Wood finishes:



Pippy Oak (PO)



Beech (B)



Odessa Golden Oak (OGO)



Code	Description	Height (mm)	Width (mm)	Depth (mm)
CTS	Square coffee table	630	460	460
CTCL	Circular coffee table - low height	460	Ø600	N/A
CTC	Circular coffee table	640	Ø600	N/A
CTR	Rectangular coffee table	920	460	460
CTU	Square coffee table - low height	520	600	600

Venetian chairs

Comfortable and stylish with a robust frame, suitable for most areas.



Code	Description	Height (mm)	Width (mm)	Depth (mm)	Seat height (mm)	Seat width (mm)	Seat depth (mm)
V300	Low back chair with rounded arms	864	584	610	457	483	432
V310	Medium back chair with rounded arms	1041	610	660	457	483	482
V020	Side chair with arms	820	580	550	457	495	482



V300 Low back chair with rounded arms



V310 Medium back chair with rounded arms

Features:

- Sacral gap on all chairs for easy cleaning
- Strong and durable laminate frame
- V020 has easy to grip arms to aid rising
- V300, V080 and V020 stackable 3 high.

V020 has a curved back to ensure maximum comfort when seated for long periods.



V020 Side chair with arms

Finland chairs

A modern style available in a choice of contemporary colours.



Code	Description	Height (mm)	Width (mm)	Depth (mm)	Seat height (mm)	Seat width (mm)	Seat depth (mm)
F833	Stacking chair	850	550	530	483	500	570
F834	Stacking chair with curved arms						
F835	Stacking chair with straight arms						



F833 Stacking chair



F835 Stacking chair with straight arms

Features:

- Sacral gap on all chairs for easy cleaning
- Strong and durable laminate frame
- Stackable 6 high
- Cushioned seat and back rest.



F834 Stacking chair with curved arms

Waiting area tub chairs

We have an extensive range of tub chairs suitable for all NHS environments. With various designs and styles, they are perfect for sprucing up the waiting area.



Chelford

A classic design with a staple, smooth, curved back feature.

Features:

- Excellent durability and comfort
- Available in a variety of anti-microbial and waterproof upholstery for ease of cleaning
- Stylish, classic design.



CHELT Tub chair



Code	Description	Height (mm)	Width (mm)	Depth (mm)	Seat height (mm)
CHELT	Tub chair	890	780	720	406

Kingsley

The modern and robust design helps give you a stylish, yet practical easy chair.

Features:

- Excellent durability and comfort
- Available in a variety of anti-microbial and waterproof upholstery for ease of cleaning
- Contemporary design
- Seat height options 18" or 20".



KINL Low back chair



Code	Description	Height (mm)	Width (mm)	Depth (mm)	Seat height (mm)
KINL	Low back chair	838	610	711	457 - 508

Balero

Comfort and support with a clean and practical design.

Features:

- Sacral gap allows easy cleaning
- All wooden frames are coated in antibacterial lacquer and anti-microbial fabrics are available
- Seat height standard 19"
- Seat height optional 16", 17", 18" or 20".



BL105 Low back chair



Code	Description	Height (mm)	Width (mm)	Depth (mm)	Seat width (mm)	Seat depth (mm)
BL105	Low back chair	940	610	635	483	483

Seat height options: 16", 17", 18", 19", 20" - (19" as standard)

Teramo

Teramo chairs have been cleverly created to look good in any location.

Features:

- Optional pressure reducing cushion
- Excellent durability and comfort
- Available in a variety of anti-microbial and waterproof upholstery for ease of cleaning
- Contemporary twist on a classic design
- Seat height options 18" or 20".



TERL Low back chair



Code	Description	Height (mm)	Width (mm)	Depth (mm)	Seat height (mm)
TERL	Low back chair	840	670	640	457 - 508

Napoli

A stylish contemporary chair. Fully upholstered for superior comfort.

Features:

- Integrated seat cushion
- Sturdy arms aiding push up from chair
- Available in a variety of anti-microbial and waterproof upholstery for ease of cleaning
- Standard seat height 18".



NAP Tub chair



Code	Description	Height (mm)	Width (mm)	Depth (mm)	Seat height (mm)
NAP	Tub chair	860	590	590	457

Modern style occasional chair, with smooth arms, clean lines and unique shape.

Features:

- Integrated seat cushion
- Soft scooped arms aiding push up from chair
- Available in a variety of anti-microbial and waterproof upholstery for ease of cleaning
- Standard seat height 18".



SHRT Tub chair



Code	Description	Height (mm)	Width (mm)	Depth (mm)	Seat height (mm)
SHRT	Tub chair	780	590	500	457

Seating that makes you feel close to home.

Our range of seating has been carefully selected to withstand day-to-day demands. Every piece of furniture in this range is constructed using strong, durable materials with many featuring built-in extras.

- 40 Belton
- 41 Forli
- 42 Kingsley
- 43 Harvard
- 44 Chelford
- 45 Modica
- 46 Teramo
- 47 Milan



Belton

Stylish looks and modern practicality, in a high quality leather effect finish.




Code	Description	Height (mm)	Width (mm)	Depth (mm)	Seat height (mm)	Seat width (mm)	Seat depth (mm)
BEL	Armchair	940	910	970	480	560	530
BEL2S	2 seater sofa	940	1470	970	480	1120	530
BEL3S	3 seater sofa	940	2010	970	480	1650	530



BEL2S 2 seater sofa



FOR2S Compact 2 seater



Code	Description	Height (mm)	Width (mm)	Depth (mm)	Seat height (mm)	Seat depth (mm)	Seat width (mm)
FOR	Chair	820	730	710			510
FOR2S	Compact 2 seater	820	1280	710	457	490	1040
FOR3S	Compact 3 seater	820	1580	710			1340

Forli

Affordable, contemporary and comfortable, compact seating suitable for any living space.



BEL Armchair



FOR Chair

Kingsley

Sophisticated and simple shapes help create a stylish look.



Code	Description	Height (mm)	Width (mm)	Depth (mm)	Seat height (mm)	Seat width (mm)	Seat depth (mm)
KINL2S	Low back 2 seater	838	1219	711		455	
KINH	High back chair	1041	610	711	457 - 508	455	485
KINL	Low back chair	838	610	711		1020	



KINL2S Low back 2 seater




KINH High back chair



KINL Low back chair

Harvard

Timeless elegance with a practical and homely feel.



Code	Description	Height (mm)	Width (mm)	Depth (mm)	Seat height (mm)	Seat width (mm)	Seat depth (mm)
HAC	Armchair	910	940	890	432	483	559
HA2S	2 seater sofa	910	1410	890	432	940	559
HA3S	3 seater sofa	910	1890	890	432	1423	559



HA2S 2 seater sofa



HAC Armchair

Chelford

Flexible options come in a contemporary shape.



Code	Description	Height (mm)	Width (mm)	Depth (mm)	Seat height (mm)	Seat width (mm)	Seat depth (mm)
CHELL	Low back chair	780	720	890	457 - 508	520	432
CHELT	Tub chair	890	780	720	432	520	500
CHEL	High back chair	1060	780	720	457 - 508	520	500



CHELT Tub chair



CHEL High back chair



CHELL Low back chair

Modica

Comfortable and traditional with a robust frame.



Code	Description	Height (mm)	Width (mm)	Depth (mm)	Seat height (inches)	Seat width (mm)	Seat depth (mm)
MODW	High back wing chair	1160	830	700		450	
MODL	Low back chair	840	760	640	457 - 508	450	410
MODW2S	High back wing 2 seater	1160	1440	700		1030	



MODW2S High back wing 2 seater



MODW High back wing chair



Teramo

Teramo chairs have been cleverly created to look good in any location.



Code	Description	Height (mm)	Width (mm)	Depth (mm)	Seat height (inches)	Seat width (mm)	Seat depth (mm)
TERW	High back wing chair	1195	720	700		450	
TERL	Low back chair	840	570	640	457 - 508	450	410
TERW2S	2 Seater	1195	1330	700		1030	



TER2WS 2 Seater



TERW High back wing chair



TERL Low back chair

Milan

Eye-catching, attractive seating with flexible and practical adaptability.



Code	Description	Height (mm)	Width (mm)	Depth (mm)	Seat height (mm)	Seat width (mm)	Seat depth (mm)
MIL2SB	2 seater sofa bed	800	1562	1829	419	1500	495



MIL2SB 2 seater sofa bed (open)



MIL2SB 2 seater sofa bed

Dining furniture that makes you feel close to home.


The items in the comprehensive dining range come in various shapes to ensure whatever the need, it can be met. Our dining chairs offer flexible options depending upon the different requirements of the end user.

- 50 Dining tables
- 51 Ridley table
- 52 Chelford chair
- 53 Finland chair



Dining tables

A variety of stylish shapes and designs to cover all dining requirements, in a smooth, heat resistant laminate finish.



Code	Description	Height (mm)	Width (mm)	Depth (mm)
DTC	Circular dining table	760	1070	N/A
DTS	Square dining table		910	910
DT2S	Rectangular dining table (2 seater)		914	610
DTR	Rectangular dining table		1520	910



DTC Circular dining table

Wood finishes:



Standard - Pippy Oak (PO)

Optional - Odessa Golden Oak (OGO), Beech (B)



DTS Square dining table



DT2S Rectangular dining table (2 seater)



DTR Rectangular dining table

Ridley table

Contemporary chrome pedestal based dining table, suitable for both dining and bistro environments.



Code	Description	Height (mm)	Width (mm)	Depth (mm)
RIDDTC	Circular pedestal dining table	755	N/A	Ø1070
RIDADDTC	Adjustable circular pedestal dining table			

Wood finishes:



Pippy Oak (PO)




Odessa Golden Oak (OGO)



Fireside Select Walnut (FSW)



Beech (B)



Chrome pedestal base




18mm MDF top and sprayed bullnose edge




RIDDTC Circular pedestal dining table

Chelford chair

Contemporary and versatile with a variety of options.



Code	Description	Height (mm)	Width (mm)	Depth (mm)	Seat height (mm)	Seat width (mm)	Seat depth (mm)
CHE	Dining chair	889	495	470	457	457	432
 CHEWA	Dining chair with arms	889	495	470	457	457	432



 CHEWA Dining chair with arms

Features:

- Modern and versatile dining chair
- Available in a choice of antimicrobial and waterproof finishes
- Available with and without arms.

CHE Dining chair

Finland chair

A modern style available in a choice of contemporary colours.




Code	Description	Height (mm)	Width (mm)	Depth (mm)	Seat height (mm)	Seat width (mm)	Seat depth (mm)
F833	Stacking chair	850	550	530	483	500	570
F834	Stacking chair with curved arms	850	550	530	483	500	570
 F835	Stacking chair with straight arms	850	550	530	483	500	570



F833 Stacking chair



 F835 Stacking chair with straight arms

Features:

- Sacral gap on all chairs for easy cleaning
 - Strong and durable laminate frame
- N.B. Recommended stacking six chairs maximum.

F834 Stacking chair with curved arms

Specialist seating that makes you feel close to home.

With over 50 years' experience, we ensure that all care home reclining chairs we offer are not only comfortable but also designed and built to the highest quality.

- 56 Recliners
- 57 Siena
- 58-59 Horizon



Renray Recliner



RR2 Manual recliner

RR1 Single three way electric rise recliner



An extremely comfortable recliner available in a manual or electric version.

Features:

- Braked castors (BC)
- Pressure reducing cushion (P)
- Back height 762mm
- Features a waterfall style back.



The waterfall back available on the Renray recliner is fibre filled and can be adjusted to suit the user's needs. This allows for excellent support in any position, whether seated or reclined.



Code	Description	Seat height (mm)	Seat width (mm)	Seat depth (mm)	Max. weight
RR1	Single three way electric rise recliner	457	508	508	120kg/19st
RR2	Single three way manual rise recliner				
RR3	Manual recliner		559		159kg/25st

Optional	Code	Description
	V	Pressure cushion - visco
	PG	Pressure cushion - gel
	BC	Braked Castors

Siena



RSM2 Manual recliner



RSP1T25 Electric rise and decline tilt in space, 1 motor 2 way

A modern, sleek design recliner with clean lines and tilt-in-space function.

Features:

- Assists with sit-to-stand position
- Available in manual and electric
- 25 stone versions have tilt-in-space function
- Designed to be easy to clean.



Code	Description	Height (mm)	Width (mm)	Depth (mm)	Seat height (mm)	Seat width (mm)	Seat depth (mm)	Max. weight
RSP1	Electric rise and decline, 1 motor 3 way	1900	700	800	455	535	510	120kg/19st
RSP2	Electric rise and decline, 2 motor 3 way							120kg/19st
RSP1T25	Electric rise and decline tilt in space, 1 motor 2 way							159kg/25st
RSM2	Manual recliner, 3 way							120kg/19st
RSMT25	Manual recliner tilt in space, 2 way							159kg/25st
RSP325	Bariatric, riser recliner, 3 motor 3 way							159kg/25st
RSP340	Bariatric, riser recliner, 3 motor 3 way	255kg/40st						

Options	Code	Description
	BC	Braked castors
	V	Pressure reducing visco cushion
	PG	Pressure cushion - gel

Horizon

Features:

- Easy to grip to assist getting in, up and out of chair
- Fixed back for controlled support
- Removable duvet seat pad for cleaning or updating
- Clean-lined zipped back cover with easy access to motor
- Integrated pressure relief seat pad options
- Standard castors for easy manoeuvrability
- Accessories - lumbar cushions, 10 day battery backup, USB port.



Code	Description	Seat height (mm)	Seat width (mm)	Seat depth (mm)	Max. weight	Options
RHP1	1 motor rise recline	457	508	508	120kg/19st	BC Braked castors
RHPT25	1 motor tilt in space				159kg/25st	V Pressure reducing visco cushion
RHM6	Manual recline only				120kg/19st	PG Pressure cushion - gel
RHMT25	Manual tilt in space recline				159kg/25st	
RHP340	3 motor rise recline	457	737	508 - 610	255kg/40st	

Soft flowing, clean-lined, unfussy, modern design. The defined top-stitch detailing enhances the contoured styling and smart proportions of the design. Simply designed with maximum appeal and uncompromised comfort, Horizon is easy on the eye, easy to accommodate, easy to blend in.



Beds that make you feel close to home.

With over 50 years' experience, we can guarantee that all our beds are designed and built with you in mind. From comfort to assembly, we pride ourselves on offering top quality products.

- 62 Elite
- 63 Prestige bariatric
- 64 Carer
- 65 Corus
- 66-68 Mattresses
 - 66 - Corus mattress
 - 67 - Nautilus mattress
 - 67 - Bariatric mattress
 - 68 - Jet mattress
- 69 Seat cushions
 - 69 - Corus foam
 - 69 - Noodles technology
 - 69 - Aurora active
- 70-73 PAC
 - 70-71 - Noodles mattress
 - 72 - Dormir
 - 73 - Diamond
- 74-77 Beds
 - 74 - Abbey divan bed
 - 75 - Denver box bed
 - 76 - Jetbed
 - 77 - Stretford divan bed



Elite

The Elite profiling bed, complete with full length wooden side rails. Standard height range 40-80cm, suitable for residents and carers to use safely. Side rail lock-down facility enabled with a key to safely secure side rails when specified. A padded headboard is available in a choice of fabrics. This 4 section bed arrives disassembled in transport mode and can be assembled by one person. Mattress NOT included.



Features:

- Overall length 2200mm
- Height adjustment 275 - 675mm
- Lockable handset controls all electrically operated functions
- Lockable castors
- Low side rails for easy access
- Side rail lock (optional)
- Full length wooden side rails
- Contrasting bumper
- Rounded corners
- Low profiling height 275mm
- Head and foot down tilt
- Curved head end
- Dementia-friendly.



Code	Description	Height (mm)	Width (mm)	Length (mm)	Max. weight
E2000LKPH	4 section bed including full length wooden side rails and key side rail lock down facility and padded headboard	275 - 675	1040	2200	180kg/28st
E2000LK	Elite 4 section bed including full length wooden side rails and key side rail lock down facility				

Prestige bariatric

Prestige bariatric low height bed to suit up to 50 stone (320kgs) person, a 4 section electric profiling bed with side rails.

Features:

- Robust and sturdy electrically profiling bariatric bed for use in the Community, Nursing and Residential Care markets
- Quick and easy to assemble, even within limited spaces
- Low platform height range of 24.5cm to 61.5cm
- Lockable handset controls all electrically operated functions
- Four braked castors for added safety
- Electrically operated rising backrest and kneebreak functions
- Infinitely variable tilt positioning.



Code	Description	Height (mm)	Width (mm)	Depth (mm)	Max. weight
PTGBAR2	Bariatric low bed, 4 section bed including side rails	275 - 690	1340	2195	380kg/60st


Carer

A cost-effective solution to compliance, safety and care. Carer 4 section bed including full length wooden side rails, lock down facility. Variable height range 19.5 - 60cm. Mattress NOT included.



Features:

- IEC 60601-2-52 compliant
- 4 section electric profiling bed - easy transportation, assembly and installation by a single person
- Exceptional low height of 19.5cm - suitable for falls risk users
- Enhanced profiling system for improved patient positioning with reduced shear and friction
- Tilt head down to -10.5 degrees and reverse tilt (head up) to +10.5 degrees
- Integral bed extension of 11.5cm with mattress retainer
- Vascular and fowler leg positioning
- Battery backup/emergency lowering facility.



Code	Description	Height (mm)	Width (mm)	Depth (mm)	Max. weight
CAR1000K	4 section bed including full length wooden side rails, lock down facility	195 - 600	1040	2207	180kg/28st

Corus

Corus low profile bed has a unique low height of 19.5cm and is designed to help in the prevention of falls and injuries to clients. The simple to use variable height allows greater independent movement and transfers.



Features:

- IEC 60601-2-52 compliant
- 4 section electric profiling bed - easy transportation, assembly and installation by a single person
- Exceptional low height of just 19.5cm - suitable for fall-risk users
- Enhanced profiling system for improved patient positioning with reduced shear and friction
- Tilt head down to -10.5 degrees and reverse tilt (head up) to +10.5 degrees
- Integral bed extension of 11.5cm with mattress retainer
- Vascular and fowler leg positioning
- Battery backup/emergency lowering facility.



Code	Description	Height (mm)	Width (mm)	Depth (mm)	Max. weight
COR1000	Low height profiling bed. 3 section electrical profiling bed. 100mm kick stop braking castors, complete with transport brackets	195 - 600	937	2163	180kg/28st

Corus mattress

The Corus foam mattress is an engineered foam, pressure redistributing mattress, providing high levels of protection and comfort for your residents.



C4000 Very high risk pressure reducing foam mattress




Corus 4000 mattress

Features:

- Designed to contour with electrical profiling bed
- 3 zone castellated foam offering excellent pressure redistribution
- Reinforced sides which provides excellent stability
- Vapour permeable waterproof PU cover with concealed zip on 3 sides.



C5000 Very high risk REVO gel infused pressure reducing foam mattress



Corus 5000 mattress

Features:

- REVO gel infused visco elastic top layer insert to aid with temperature control
- Firm edges for safer egress and firm base layer foam
- Modular U-shaped base layer
- Two way turn
- Waterproof PU cover with concealed zip on 3 sides.


Code	Description	Height (mm)	Width (mm)	Depth (mm)	Max. weight
C4000	Very high risk pressure reducing foam mattress.	150	850	1970	240kg/38st
C5000	Very high risk REVO gel infused pressure reducing foam mattress.				180kg/28st

Nautilus mattress

An exceptional waterproof mattress that is hypo-allergenic, antimicrobial and breathable.



NM Mattress (single)



NMD Mattress (double)

Code	Description	Height (mm)	Width (mm)	Depth (mm)
NM	Mattress (single)	175	915	1905
NMD	Mattress (double)	200	1350	1900

Features:

- 100% transfer printed polyester, impregnated and multi-coated with micro-porous polymers and flame retardant
- Easy to clean; water-based and oil based stains can be removed with household cleaning agents.

Bariatric mattress

Bariatric pressure reducing mattress.

Code	Description	Height (mm)	Width (mm)	Depth (mm)	Max. weight
BHCMATT	Bariatric pressure reducing mattress				320kg/50st

Features:

- Maximum user weight 320kg (50 stone)
- For use with the PTGBAR2 bariatric bed.



BHCMATT Pressure reducing mattress. For use with PTGBAR bariatric bed


Jet mattress

The Jet mattress is a pressure redistributing mattress designed to be the first line of defence against tissue damage in Mental Health Environments.



Features:

- Conforms to all bed profile positions
- Castellated cut foam for better pressure redistribution
- Vapour permeable cover for enhanced micro-climate
- Fully welded seams for enhanced infection control
- Sealed mattress with no access for contraband
- CRIB 7 Fire Rated
- Designed for CE Environments
- Tough and long lasting without clinically compromising the patient.



Code	Description	Height (mm)	Width (mm)	Depth (mm)
JETMATPR	Pressure reducing Crib 7 vandal proof mattress waterproof permeable PU welded seams	416	930	1937
JETMATT5	Crib 7 vandal proof mattress waterproof permeable PU welded seams	175	900	1900

Seat cushions

Our range of foam, non-foam and active cushions provide a range of seating options to provide pressure relief and contours around the user's body to provide both support and comfort.

Corus foam

Code	Description	Height (mm)	Width (mm)	Depth (mm)	Max. weight
CC3000	High risk, foam, pressure reducing cushion with contoured front edge zoned to surface. Non slip pad on base. Vapour permeable, waterproof PU cover with concealed zip on 3 sides.	75	430	430	89kg/14st
CC5000	Very high risk foam pressure reducing cushion with contoured front edge. Visco elastic foam top surface. Non slip pad on base. Vapour permeable, waterproof PU cover with concealed zip on 3 sides.				



Noodles technology

Code	Description	Height (mm)	Width (mm)	Depth (mm)	Max. weight
NOODCO-90	Noodles technology cushion high risk non-foam cushion fully cleanable core and cover providing greater immersion and improved micro climate	75	430	430	110kg/17st



Aurora active

Code	Description	Height (mm)	Width (mm)	Depth (mm)	Max. weight
AUR046046D	Alternating air cushion with removable PU water proof/ vapour permeable cover complete and advanced "whisper quiet" pump. Manual comfort setting. For high risk residents.	75	460	460	108kg/17st



Features:

- All our cushions reduce or redistribute pressure
- Suitable for high risk users
- Designed with a removable washable cover
- Covers are made of waterproof PU material
- Concealed zip for better infection control.

Noodles mattress


The Noodles Technology mattress range is a revolutionary, fully cleanable and recyclable, zoned static mattress system, designed to address the challenges of those that are bed bound or clinically require a static mattress.



Code	Description	Height (mm)	Width (mm)	Depth (mm)	Max. weight
NOODM090200-180	Mattress very high risk non-foam system fully cleanable core and cover U-shaped base with firm edges zoned insert with heel section providing greater immersion and improved micro climate.	850	1870	150	180k/28st

Benefits:


- Fully cleanable**
 Improving infection control issues
 - Improved micro-climate**
 Reduces Micro-Climate humidity by 54% and temperature by 17% due to open Noodles structure and Dartex cover
 - Recyclable**
 Improving carbon footprint
 - Enhanced patient comfort**
 Superb all-body support, helping to reduce the risk of pressure ulcers
 - Cost effective solution**
 Longer life product, meaning reduced replacement costs
- The NOOD085197/40:25 mattress is typically used in community loan stores and acute hospital wards.



Dormir



The Dormir active mattress replacement system, provides high levels of care and comfort with an ultra-low 12.5cm profile.



Code	Description	Height (mm)	Width (mm)	Depth (mm)
 DOR090200BD	Dormir mattress DOR303 and PUMP8535 high risk low profile air mattress system with 29 x 10cm compact cells and 3cm foam base. Removeable PU top cover, CPR. Digital comfort setting.	125	880	1980

Features:

- Low profile 12.5cm, no need to use side rail height extensions
- 29 compact cell construction for best possible pressure area care
- Whisper-quiet pump ensures tranquil bedrooms
- PU cell design which provides durability and minimises friction and shear when the system is profiled.




Dormir cell technology.

Diamond


The Diamond is an active mattress replacement system, providing the highest level of care and comfort.



Code	Description	Height (mm)	Width (mm)	Depth (mm)	Max. weight
 DIA090200	Diamond standard 200cm long, 90cm wide, 12.5cm low profile, dynamic hybrid air mattress replacement system complete with removable PU water proof/vapour permeable cover, CPR facility and advanced "whisper quiet" pump with visual/audible low pressure and power failure alarms suitable high/very high risk residents or those at risk of falls	125	880	1980	240kg/28st

Features:

- 12.5cm ultra-low profile, does not need side rail height extensions
- Stacked cell technology, giving the best possible pressure area care
- Whisper-quiet pump for tranquil bedrooms.



Diamond cell technology.

Abbey divan bed

Effortlessly stylish and practical, the Abbey divan beds are the perfect solution for a peaceful night's sleep.

Features:


- Includes head and foot board
- Base only, mattress NOT included
- SWL 165kg
- Max patient weight 150kg.



BEDAB Single bed (base only), including head and foot board



ABEDDS Bed with dimple base




Code	Description	Height (mm)	Width (mm)	Depth (mm)
BEDAB	Single bed (base only), including head and foot board	880	970	2025
ABEDDS	Bed with dimple base			

Denver box bed

The Denver box bed is designed to accommodate the needs of a more challenging environment.

Features:

- 18mm thickness
- Security screws.



Code	Description	Height (mm)	Width (mm)	Depth (mm)
DENS8971CE	Box bed for challenging environment	460	971	1962

Jet bed

The Jet bed has been specifically designed to provide a furniture solution for mental health facilities to combat cleaning and security problems associated with box beds, foam beds and more traditional style beds.



JETBASE 21" Deep crib 5 foam base to be used under jet mat



JETBEDCE Jet bed challenging environment

Features:

- Profiled leg
- Head and foot panels integral
- Reduced lip
- Easier to move compared to box beds
- 18mm solid wood internal structure
- All panels laminated MDF
- All edges in floor facing lacquered
- No PVC taped edges
- Optionally sealed during construction
- No visible screws.

Benefits:

- Reduced ligature risk when bed stood on end
- Reduced risk of head and foot panels being 'kicked' off as part of the solid structure
- Easy cleaning, infection control friendly
- Enables deep cleaning and security inspections
- Solid, heavy and robust
- Hard wearing, scratch resistant and fully wipeable
- Prevent water ingress, easy cleaning, tamper proof
- More durable, no fluid ingress, tamper proof
- Less risk of 'picking', more reliable solution
- No risk of tampering, safe and secure.



Code	Description	Height (mm)	Width (mm)	Depth (mm)
JETBEDSCE	Silicone sealant for challenging environment	416	930	1937
JETBASE	21" Deep crib 5 foam base to be used under jet mat.	416	930	1937

Stretford divan bed

A traditional, practical and robust style divan bed suitable for healthcare environments.



BEDDB Double divan bed with nautilus waterproof fabric



BEDSB Single divan bed with nautilus waterproof

Features:

- Very strong and durable
- 6 sturdy wooden legs
- Crib 5 fire resistant and waterproof cover
- Optional headboard.



Code	Description	Height (mm)	Width (mm)	Depth (mm)
BEDSB	Single divan bed with nautilus waterproof fabric (Base only)	345	910	1910
BEDDB	Double divan bed with nautilus waterproof fabric (Base only)	400	1350	1900

Staff accommodation that makes you feel close to home.

Our contemporary and stylish bedroom furniture ranges create a calm and homely environment to help staff relax and unwind after a long shift.



80 Florida range
81 Albany range




Florida range



FL315 Double wardrobe

The Florida range is simple without compromising on design, creating a sophisticated and contemporary look in its simplicity.



Code	Description	Height (mm)	Width (mm)	Depth (mm)
FL116	Bedside unit with door and lock	660	445	457
FL117	3 drawer bedside unit and lock	660	445	457
FL245	30" 4 drawer chest	610	762	457
FL315	Double wardrobe	1829	762	584
FL415	Top box	400	762	584
FL722	Sideboard - 2 door, 2 drawer	840	900	500
FL761	Bookcase - floor standing	1000	900	350

Features:

- 18mm MFC drawer and door fronts with PVC edging
- 18mm MFC tops with PVC edging
- Chrome bar handles
- Lockable bedside unit - FL116, FL117
- Available in a variety of wood finishes.

Wood finishes:



Pippy Oak (PO)



Odessa Golden Oak (OGO)



FL415 Top box



FL116 Bedside unit with door and lock



FL117 3 drawer bedside unit and lock



FL245 30" 4 drawer chest



FL722 Sideboard - 2 door, 2 drawer



FL761 Bookcase - floor standing

Albany range



ALB315 Double wardrobe

The Albany range is refined, practical and stylish, a stunningly simple option for any bedroom environment.



Code	Description	Height (mm)	Width (mm)	Depth (mm)
ALB116L	Bedside unit with door and lock	650	450	432
ALB255	5 drawer chest	855	815	432
ALB315	Double wardrobe	1800	815	610

Features:

- 18mm MFC drawer and door fronts with PVC edging
- 18mm MFC tops with PVC edging
- Chrome D handles
- Inframe design
- Available in a variety of wood finishes.

Wood finishes:



Oxide Vintage Oak (OVO)



Odessa Golden Oak (OGO)



Oxide Vintage Oak and Cashmere fronts (OCO&C)



ALB116L Bedside unit with door and lock



ALB255 5 drawer chest

High performance fabrics for a ward environment

Amongst our extensive range of vinyl and impervious fabrics, we offer a number of anti-microbial options that are effective against MRSA as well as E.coli, Listeria, Salmonella, Legionella and other similar bacteria.

Benefits of vinyl fabrics include: waterproof, anti-microbial, anti-fungal, bleach cleanable, flame retardant and scratch resistant. Our fabrics are suitable to be cleaned with the following products: Antichlor, Antichlor Plus, Haz Tabs, Chlor Clean, 70% Alcohol Wipes (Caveat: Check with your IPC department which cleaning products you are allowed to use).



For sample swatches or copies of the cleaning/care instructions on your chosen fabrics please contact our Sales Office on **01606 593456**.



Cleaning Guide

Source: http://www.panaz.com/media/site-downloads/CARE_MAINTENANCE_OF_UPHOLSTERY_updated_.pdf

Provided by Panaz Fabrics


Care and maintenance of vinyl upholstery*

01. Wipe clean regularly with a microfibre cloth to remove dust particles and light soiling.
02. Remove spillages promptly with an absorbent dry cloth.
03. For general soiling, use a microfibre cloth (taking care to avoid abrasive damage to surface) and tepid soapy water (non-alkaline).
04. Rinse clean with water and dry with an absorbent cloth.
05. Bodily fluids should be removed promptly within 15 minutes with cold water and then cleaned as above.
06. Extreme soiling/infection control purposes use a 5% sodium hypochlorite bleach solution, Haz Tabs or Chlorclean.
07. Steam cleaning is an effective sanitizer as are Haz Tabs and Chlorclean but is not a good cleaner. A wipe with a cloth, preferably microfibre, and water is more effective for actual cleaning. Steam and heat from steam cleaners actually seal in oily residues from bodily fluids, etc. If steam cleaning is essential, use infrequently for emergency infection control purposes only and in accordance with the machine manufacturer's recommendations. As we are unable to comment on the variety of machines in the marketplace or on the operator's handling, care should be taken for the sake of both the upholstery fabric and the wooden frame and we are unable to guarantee these.
08. Do not machine wash or dry clean.
09. Do not use furniture polish, stain remover or solvents.
10. Remove indigo dye from denim and canvas garments immediately as they will produce permanent staining if not dealt with promptly.
11. Do not use a felt tip or ball-point pen to mark.

*Refer to fabric manufacturer's guidelines

Express delivery


Certain items are available for express delivery from stock and can be with you within 15 business days from receipt of order. For these Express Delivery items, look out for the special icon which is shown next to the relevant products.




LOOK OUT FOR THE EXPRESS DELIVERY ICON ON MANY OF OUR PRODUCTS.

Express seating products are supplied in the following vinyls and wood finish:


Upholstery vinyl




105 Indigo



198 Hyacinth




110 Pacific



837 Stone

Wood



Beech (B)

A full turnkey solution

Renray are specialists in providing a complete turnkey package, ranging from chairs and tables, patient bedside furniture plus specialised seating, beds and mattresses. The package also includes free room placement, free delivery, full project management, relevant product training and removal of waste materials.

Visit our showrooms

To help you make your decision, visit our showrooms where you can view and try out our furniture. The showrooms are located at our head office, just off junction 18 of the M6 in Winsford, Cheshire.

To arrange a visit please contact our Sales Office on **01606 593456**.

Order delivery

Delivery is generally between 4 to 6 weeks after receipt of your order. Please enquire if a quicker delivery is required (see page 83 for Express Delivery details).

All orders are delivered by our own fleet of lorries and wagons and installed by our experienced, DBS checked uniformed staff.

Quotations and samples

For a competitive quotation, please contact our Sales Office on 01606 593456. One of our nationally based Area Sales Managers will be happy to arrange a meeting at a location of your choice to discuss your needs. Furniture samples can be made available to you and delivered for a two week trial period.

After sales service

If you have any further needs after your order has been delivered, or any comments and suggestions, please contact us. We value the ongoing relationship we have with our customers.



Renray credentials

When you work with Renray, you're in very good company and here's why...

Certifications

When ordering any Renray product, you have the peace of mind of knowing that we have been audited and accredited and have the necessary standards and business controls in place. These are the professional business standards when dealing with at risk and vulnerable clients.



ISO 9001 - Quality Management System



ISO 14001 - Environmental Management System



ISO 18001 - Health and Safety Management System



FSC - FSC® - Forestry Stewardship Council®
Our policy is to provide and maintain working practices that have the least detrimental effect on the global environment. We operate a lean manufacturing process whereby all our products are made to order by our highly skilled workforce. Look out for our FSC®-certified products.



CHAS - The Contractors Health and Safety Assessment Scheme



BCFA - British Contract Furnishing Association

Professionally made

With over 50 years' experience behind us, our craftsmen consistently produce quality furniture and interiors from our own factories, based in the UK and Europe. Together, these total over 100,000 sq.ft of production space.

Renray products are independently tested and manufactured in accordance with recognised British standards.



These nationally recognised standards include:

- BS 5852 Part 2, Crib 5 Flame Retardancy
- BS EN 16139:2013 Strength, Durability & Safety. Non-domestic Seating
- BS EN 15372:2008 Strength, Durability & Seating. Non-domestic tables
- BS EN 16121:2013 Non-domestic storage. Safety, Strength, Durability & Stability
- BS EN 7176:2007 Resistance to Ignition
- BS 5867-2 type B & C Fabric for curtains, drapes & blinds
- BS 7177:2008 Resistance against ignition of mattresses, mattress pads, divans and bed bases.

Carefully installed

Renray have an in-house logistics team who co-ordinate and manage a weekly delivery service to all areas of the UK and Ireland, from our 38 swap body fleet.

Products are delivered by our own team of experienced DBS checked, uniformed drivers.

Free room placement, unpacking and packaging removal.

Highly trained staff with many years experience in the careful treatment of installing turnkey solutions.

Contact

Your next step...

For further information, including product availability, the latest colour options and prices, and to order any of the products in this brochure; or if you have any queries about other Renray products and services, contact our Sales Office on:



01606 593456



sales@renrayhealthcare.com



01606 861354

Our experienced team will be pleased to discuss your requirements. Also, our nationally based Area Sales Managers are always available to visit you for a no obligation discussion and to help advise you.



<https://www.facebook.com/renrayhealthcare/>



<https://twitter.com/RenrayHC>



<https://www.linkedin.com/company/11013192/>



Road Five, Winsford Industrial Estate, Winsford, Cheshire CW7 3RB
T: 01606 593456 F: 01606 861354 E: info@renrayhealthcare.com

renrayhealthcare.com